

Product datasheet for **MP217870**

Ubr5 Mouse qPCR Primer Pair (NM_001081359)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	70790
Forward Sequence:	CGTCCTTCAGGACTGGAATGCT
Reverse Sequence:	CTCTGCTGGTTGAGGTAGACCT
ACCN:	NM_001081359 , NM_001081359.1 , NM_001081359.2 , NM_001081359.3 , BC034075 , BC037687 , BC047348 , BC049162 , BC049224 , BC051104 , BC057458 , BC057923 , BC064751
UniProt ID:	Q80TP3
Synonyms:	Edd; Edd1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)