

Product datasheet for **MP217993**

Ugt1a9 Mouse qPCR Primer Pair (NM_201644)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	394434
Forward Sequence:	TTGGTGGGATCAACTGCCTCCA
Reverse Sequence:	CGGAATCTCTGAGACCATGGATC
ACCN:	<u>NM_201644</u> , <u>NM_201644.1</u> , <u>NM_201644.2</u> , <u>BC138699</u> , <u>BC019434</u> , <u>BC146021</u>
UniProt ID:	<u>Q62452</u>
Synonyms:	mUGTBr/P; UDPGT; UDPGT 1-9; UGT1-09; UGT1.9; Ugt1a12; UGTP4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)