

Product datasheet for **MP218075**

Sun1 Mouse qPCR Primer Pair (NM_024451)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	77053
Forward Sequence:	ATCCGAGGCTATTGTGTCTGCC
Reverse Sequence:	CTCCAGAGCAAAGTCCACCATC
ACCN:	BC047928 , BC048156 , NM_024451 , NM_024451.1 , NM_024451.2 , BC030330 , BQ551727 , BY036585
UniProt ID:	Q9D666
Synonyms:	4632417G13Rik; 5730434D03Rik; mKIAA0810; Unc84a
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)