

## Product datasheet for **MP218129**

### Utp14a Mouse qPCR Primer Pair (NM\_028276)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	72554
Forward Sequence:	TGCTGAGCCACTGGGTAACCAA
Reverse Sequence:	AGACCAGGTCAACAGCCTTGC
ACCN:	<a href="#">BC082590</a> , <a href="#">NM_028276</a> , <a href="#">NM_028276.1</a> , <a href="#">BC062165</a>
UniProt ID:	<a href="#">Q640M1</a>
Synonyms:	2700066J21Rik; JsdX; mKIAA0266; NY-CO-16; SDCCAG16; T25628
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)