

Product datasheet for **MP218136**

Uvr^{ag} Mouse qPCR Primer Pair (NM_178635)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	78610
Forward Sequence:	CAAGCTGACAGAAAAGGAGCGAG
Reverse Sequence:	GGAAGAGTTTGCCTCAAGTCTGG
ACCN:	<u>BC034176</u> , <u>NM_178635</u> , <u>NM_178635.1</u> , <u>NM_178635.2</u> , <u>NM_178635.3</u> , <u>BB530952</u> , <u>BC034072</u> , <u>BY024873</u>
UniProt ID:	<u>Q8K245</u>
Synonyms:	9530039D02Rik; AI648770; BB124205; Uvr ^{ag} 1; Uvr ^{ag} l
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)