

Product datasheet for **MP218177**

Vangl2 Mouse qPCR Primer Pair (NM_033509)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	93840
Forward Sequence:	GCATCCTTCAGCACCTGGAGTT
Reverse Sequence:	TCGCTCACCAAGGTCCACTGTT
ACCN:	<u>BC052195</u> , <u>NM_033509</u> , <u>NM_033509.1</u> , <u>NM_033509.2</u> , <u>NM_033509.3</u> , <u>NM_033509.4</u> , <u>BC075636</u> , <u>BY320814</u>
UniProt ID:	<u>Q91ZD4</u>
Synonyms:	C530001F03Rik; loop-tail; Lootl; Lp; Lpp1; Ltap; ska17; ska; stbm; strabismus; Vangl12
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)