

Product datasheet for **MP218181**

Vegfa Mouse qPCR Primer Pair (NM_001025250)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22339
Forward Sequence:	CTGCTGTAAACGATGAAGCCCTG
Reverse Sequence:	GCTGTAGGAAGCTCATCTCTCC
ACCN:	<u>BC061468</u> , <u>NM_001025250</u> , <u>NM_001025250.1</u> , <u>NM_001025250.2</u> , <u>NM_001025250.3</u> , <u>BC022642</u> , <u>BQ554097</u> , <u>BQ832724</u> , <u>BU920197</u> , <u>BY001490</u>
UniProt ID:	<u>Q00731</u>
Synonyms:	V; Veg; Vegf; VEGF12; VEGF16; VEGF18; Vpf
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)