

Product datasheet for **MP218203**

Vcan Mouse qPCR Primer Pair (NM_001081249)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	13003
Forward Sequence:	GGACCAAGTTCCACCCTGACAT
Reverse Sequence:	CTTCACTGCAAGGTTCTCTTCT
ACCN:	<u>NM_001081249</u> , <u>NM_001081249.1</u> , <u>BB009616</u> , <u>BB104197</u> , <u>BC021652</u> , <u>BC096495</u> , <u>BM251152</u> , <u>BQ828249</u>
Synonyms:	5430420N07Rik; 9430051N09; Cspg2; DPEAAE; hdf; NG2; PG-M; PG-M(V0); PG-M(V1)
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)