

Product datasheet for **MP218225**

Vcl Mouse qPCR Primer Pair (NM_009502)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22330
Forward Sequence:	CCTATCAAGCTGTTGGCAGTAGC
Reverse Sequence:	TGTGGCTCCAAGCCTTCCTGAA
ACCN:	<u>BC008520</u> , <u>NM_009502</u> , <u>NM_009502.1</u> , <u>NM_009502.2</u> , <u>NM_009502.3</u> , <u>NM_009502.4</u> , <u>BC008554</u> , <u>NM_009502.5</u>
UniProt ID:	<u>Q64727</u>
Synonyms:	9430097D22; AA571387; AI462105; AW545629
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)