

Product datasheet for **MP218417**

Vsig2 Mouse qPCR Primer Pair (NM_020518)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57276
Forward Sequence:	GGCGACAACCTTGCCTTGGAGT
Reverse Sequence:	CCGATCTGCTTTAGAGCCAGTG
ACCN:	<u>BC099457</u> , <u>NM_020518</u> , <u>NM_020518.1</u> , <u>NM_020518.2</u> , <u>BB869209</u>
UniProt ID:	<u>Q9Z109</u>
Synonyms:	1190004B15Rik; 2210413P10Rik; AU040484; Ctm; CTX
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)