

## Product datasheet for **MP218448**

### Wdfy3 Mouse qPCR Primer Pair (NM\_172882)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	72145
Forward Sequence:	TGGAGCAGACAGATTCCGTCCT
Reverse Sequence:	GGAGGAAAGACTGGAGAACCAG
ACCN:	<a href="#">BC058274</a> , <a href="#">NM_172882</a> , <a href="#">BC014824</a> , <a href="#">BC034123</a> , <a href="#">BC049866</a>
Synonyms:	2610509D04Rik; ALFY; AW319683; B930017C24; Bchs; BWF1; D5Ert66e; Ggtb3; mKIAA0993; ZFYVE25
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)