

Product datasheet for **MP218470**

Wdr33 Mouse qPCR Primer Pair (NM_028866)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	74320
Forward Sequence:	GGTGCTGAATGACAGGAAAGAGG
Reverse Sequence:	GATCAGCCTGTTTCTGTGCAAGC
ACCN:	<u>NM_028866</u> , <u>NM_028866.1</u> , <u>NM_028866.2</u> , <u>NM_028866.3</u> , <u>BC003939</u> , <u>BC010638</u> , <u>BC011156</u> , <u>BC019463</u> , <u>BC019526</u> , <u>BC030076</u> , <u>BC064745</u>
UniProt ID:	<u>Q8K4P0</u>
Synonyms:	1110001N06Rik; 2310011G05Rik; 2810021O11Rik; 8430413N20Rik; WDC146
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)