

Product datasheet for **MP218477**

Dcaf12 Mouse qPCR Primer Pair (NM_026893)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	68970
Forward Sequence:	GCCATCGTACAACGTCAAGTCTG
Reverse Sequence:	CTCGGATGTCATAGAACAGCAGG
ACCN:	<u>BC023852</u> , <u>NM_026893</u> , <u>NM_026893.1</u> , <u>NM_026893.2</u> , <u>NM_026893.3</u> , <u>BY338453</u>
UniProt ID:	<u>Q8BGZ3</u>
Synonyms:	1500001L20Rik; 5830424K06Rik; AA420338; AI851081; Wdr40a
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)