

Product datasheet for **MP218604**

Xlr4b Mouse qPCR Primer Pair (NM_021365)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27083
Forward Sequence:	CAGGCTATTAAGTGTTCGCCG
Reverse Sequence:	GAACACCATCCTTTGCTTCCTGC
ACCN:	<u>BC025576</u> , <u>NM_021365</u> , <u>NM_021365.1</u> , <u>NM_021365.2</u> , <u>NM_021365.3</u> , <u>BC145482</u> , <u>BC034848</u> , <u>BC138975</u> , <u>BC138979</u> , <u>BC145481</u> , <u>BC145483</u> , <u>BE627633</u> , <u>BU755708</u>
Synonyms:	Xlr4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

