

## Product datasheet for **MP218680**

### Zbtb25 Mouse qPCR Primer Pair (NM\_028356)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	109929
Forward Sequence:	ATCCTCATCAGAGGTGACTGGC
Reverse Sequence:	GGTATGGAGGTGTTGTCTCAGG
ACCN:	<a href="#">BC068249</a> , <a href="#">NM_028356</a> , <a href="#">NM_028356.1</a> , <a href="#">NM_028356.2</a>
Synonyms:	2810462M08Rik; 2900064P18Rik; AI842128; Kup; Zfp-50; Zfp50
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)