

Product datasheet for **MP218701**

Zbtb8os Mouse qPCR Primer Pair (NM_025970)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67106
Forward Sequence:	GGAGACTCTGGAGGAAGCAT
Reverse Sequence:	AAGTGAACAGCAGACTGCAG
ACCN:	BC094629 , BC125606 , BC125608 , NM_025970 , NM_025970.1 , NM_025970.2 , NM_025970.3 , BC026442 , BC037657 , BC046825 , NM_025970.4
UniProt ID:	Q505B7
Synonyms:	2010001H09Rik; 2310028N13Rik; AI851226; Arch; Archm
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)