

Product datasheet for **MP218714**

Zc3h14 Mouse qPCR Primer Pair (NM_029334)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	75553
Forward Sequence:	AACAGGCTGCAACTTCATCCGC
Reverse Sequence:	CTTCTACGCTGCTAAGCACTGG
ACCN:	BC024824 , NM_029334 , NM_029334.1 , NM_029334.2 , BB610803 , BC021797 , BC024856 , BC026610 , BC034760 , BU936221
UniProt ID:	Q8BJ05
Synonyms:	1010001P15Rik; 1700016A15Rik; 2700069A02Rik; AU014748
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)