

Product datasheet for **MP218727**

Zeb1 Mouse qPCR Primer Pair (NM_011546)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	21417
Forward Sequence:	ATTCAGCTACTGTGAGCCCTGC
Reverse Sequence:	CATTCTGGTCCTCCACAGTGGA
ACCN:	BC139768 , BC139769 , NM_011546 , NM_011546.1 , NM_011546.2 , NM_011546.3 , BC090987 , BE852657
UniProt ID:	Q64318
Synonyms:	3110032K11Rik; AREB6; BZP; MEB1; Nil2; TCF-8; Tcf8; Tcf18; Tw; ZEB; Zfh1a; Zfx1a; Zfx1ha
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)