

Product datasheet for **MP218859**

Zfp423 Mouse qPCR Primer Pair (NM_033327)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	94187
Forward Sequence:	AGGACTTCGAGTCTCTGGCAGA
Reverse Sequence:	TGCGTAGGTGACGCAACATCCT
ACCN:	NM_033327 , NM_033327.1 , NM_033327.2 , BC139028 , BB788751 , BC059234 , BC079586 , BC139030 , BF021428 , BQ180174 , BU519630
UniProt ID:	Q80TS5
Synonyms:	ataxia1; Ebfaz; mKIAA0760; nur12; Roaz; Zfp104; Znf423
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)