

Product datasheet for **MP218969**

Zfp317 Mouse qPCR Primer Pair (NM_172918)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	244713
Forward Sequence:	ATGCTCCAGAGCCTGATACACC
Reverse Sequence:	TCTCTGGGAAGAATCTAGCAGGG
ACCN:	<u>BC059211, NM_172918, NM_172918.1, NM_172918.2, NM_172918.3, NM_172918.4</u>
Synonyms:	4932416G03; D230022C05Rik; KRAB9; mszf40; Zfp67; Zfp75
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)