

Product datasheet for **MP219003**

Zfp91 Mouse qPCR Primer Pair (NM_053009)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	109910
Forward Sequence:	AGCTCTAGGTCCAAGACTGGCT
Reverse Sequence:	CGGGAGACTGATGTTCTTCTGG
ACCN:	NM_053009 , NM_053009.1 , NM_053009.2 , NM_053009.3 , BC057323 , BC064766 , BC083000 , BC156276 , BC157098 , BC168394 , BE947731
UniProt ID:	Q62511
Synonyms:	9130014I08Rik; A530054C17Rik; AL024263; AW545902; Pz; Pzf; Zfp-91
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)