

Product datasheet for **MP219303**

Psd3 Mouse qPCR Primer Pair (NM_177698)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	234353
Forward Sequence:	CTGCCAACTCTGTCTACACCAG
Reverse Sequence:	CCATTGCTGAGACTTTCTGTGGT
ACCN:	BC003498 , NM_177698 , NM_177698.4 , BC042208 , BC052365 , BC131913 , BC145641 , BM898780
UniProt ID:	Q2PFD7
Synonyms:	4931420C21Rik; AI661273; BC003498; D430018P08; EFA6D
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)