

Product datasheet for **MP219337**

Magi1 Mouse qPCR Primer Pair (NM_001029850)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	14924
Forward Sequence:	CTGACGCACAATCAAGTCGTGG
Reverse Sequence:	GGGCTCTTCTTGGAAGTGGTA
ACCN:	<u>NM_001029850</u> , <u>NM_001029850.1</u> , <u>NM_001029850.2</u> , <u>NM_001029850.3</u> , <u>NM_001029850.4</u> , <u>BB653276</u> , <u>BC006889</u> , <u>BC060273</u> , <u>BC095943</u> , <u>BC150820</u> , <u>BU702337</u>
UniProt ID:	<u>Q6RHR9</u>
Synonyms:	AIP3; Baiap1; BAP1; Gukmi1; Magi-1; MAGI1c; mKIAA4129; TNRC19; WWP3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)