

Product datasheet for **MP219338**

Osbp19 Mouse qPCR Primer Pair (NM_133885)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	100273
Forward Sequence:	GCACCTCTTATCACAAGTCAGGC
Reverse Sequence:	GCTCACGAACAGGTCTGGATGT
ACCN:	<u>BC021507</u> , <u>BC026927</u> , <u>BC089163</u> , <u>NM_133885</u> , <u>NM_133885.1</u> , <u>NM_133885.2</u> , <u>NM_133885.3</u> , <u>NM_133885.4</u> , <u>BC002157</u> , <u>BC023759</u> , <u>BC119493</u>
UniProt ID:	<u>A2A8Z1</u>
Synonyms:	2600011I06Rik; AU015843; Orp-9
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)