

Product datasheet for **MP219789**

Mapk1 Mouse qPCR Primer Pair (NM_001038663)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26413
Forward Sequence:	TCAAGCCTCCAACCTCCTGCT
Reverse Sequence:	AGCTCTGTACCAACGTGTGGCT
ACCN:	BC058258 , NM_001038663 , NM_001038663.1 , BB648734 , BC006708
UniProt ID:	P63085
Synonyms:	9030612K14Rik; AA407128; AU018647; C78273; ERK; Erk2; MAPK2; p41mapk; p42mapk; Prkm1; PRKM2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)