

Product datasheet for **MP219865**

Cdv3 Mouse qPCR Primer Pair (NM_175833)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	321022
Forward Sequence:	TGCAGTCCACTGCCAAGCATGT
Reverse Sequence:	GCACAGCATACTGGTTGTCTAGC
ACCN:	BC052770 , BC062156 , BC096479 , NM_175833 , NM_175833.1 , NM_175833.2 , BC062156.1 , BC030486 , BC031362 , BC064103
UniProt ID:	Q4VAA2
Synonyms:	2510010F10Rik; C79446; C230084J24Rik; TPP36
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)