

Product datasheet for **MP219899**

Uba1 Mouse qPCR Primer Pair (NM_001136085)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22201
Forward Sequence:	TCGTGGTCACAGACATGGACAC
Reverse Sequence:	TGGTGGCTTGTCACCTGGATGT
ACCN:	BC058630 , NM_001136085 , NM_001136085.1 , NM_001136085.2 , BB858760 , BC138200 , BC145984 , BY004827 , BY316912
UniProt ID:	Q02053
Synonyms:	A1; A1S9; Sb; Sbx; Ube-1; Ube1x
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)