

Product datasheet for **MP220107**

Eftud2 Mouse qPCR Primer Pair (NM_001109995)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20624
Forward Sequence:	GACGAACCTGGCATCCACCTCA
Reverse Sequence:	CCTTTGGAGATGGAATGTGCTGC
ACCN:	<u>BC052674</u> , <u>BC054778</u> , <u>NM_001109995</u> , <u>NM_001109995.1</u> , <u>BC012636</u> , <u>BC036344</u> , <u>BC119549</u> , <u>BQ886207</u> , <u>BY074102</u>
UniProt ID:	<u>O08810</u>
Synonyms:	116kDa; Snrp116; U5-116kD
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)