

## Product datasheet for **MP220116**

### Zbtb33 Mouse qPCR Primer Pair (NM\_001079513)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56805
Forward Sequence:	CGAAGGCTTGCTCAACTGGTTGC
Reverse Sequence:	ACGCTTGTTGGTGGCTCACTG
ACCN:	<u>NM_001079513, NM_001079513.1, BC138647, BC025136, BC027282, BC138648, BY200826</u>
UniProt ID:	<u>Q8BN78</u>
Synonyms:	AW260253; E130014G12Rik; Kaiso; Znf-kaiso
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)