

Product datasheet for **MP220134**

Abl1 Mouse qPCR Primer Pair (NM_001112703)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11350
Forward Sequence:	TGAGCAGAAAGATGCGCCTGAC
Reverse Sequence:	CGCTCATCTTCATTTAGGCTGCC
ACCN:	BC059260 , BC103770 , NM_001112703 , NM_001112703.1 , NM_001112703.2 , BC013801 , BC035521 , BM949741 , BQ769074 , BU053120 , BY152805
UniProt ID:	P00520
Synonyms:	Abl; A1325092; c-Abl; E430008G22Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)