

Product datasheet for **MP220159**

Hras Mouse qPCR Primer Pair (NM_001130443)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	15461
Forward Sequence:	TCGCACTGTTGAGTCTCGGCAG
Reverse Sequence:	TATGCTGCCGAATCTCACGGAC
ACCN:	<u>BC011083</u> , <u>NM_001130443</u> , <u>NM_001130443.1</u> , <u>BC011083.1</u> , <u>BC050837</u> , <u>BC061885</u> , <u>BI731283</u> , <u>BQ963436</u> , <u>BU924724</u>
UniProt ID:	<u>Q61411</u>
Synonyms:	c-H-ras; c-Ha-ras; c-rasHa; H-ras; Ha-ras; Harvey-ras; Hras-1; Hras1; Kras2; ras
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)