

Product datasheet for **MP220240**

Spn Mouse qPCR Primer Pair (NM_001037810)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20737
Forward Sequence:	CATCAGCACTGCTCCAGAACC
Reverse Sequence:	TTCCGTCTGTCACAGGATTGGC
ACCN:	BC100738 , BC100739 , BC100740 , BC100741 , NM_001037810 , NM_001037810.1 , NM_001037810.2 , BB850887 , BX633083
UniProt ID:	P15702
Synonyms:	A630014B01Rik; Cd43; Galgp; Ly-48; Ly48
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)