

Product datasheet for **MP220254**

Ndor1 Mouse qPCR Primer Pair (NM_001082476)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	78797
Forward Sequence:	CCAGCAGTTCTGCCAAGTGTG
Reverse Sequence:	TACTGGGACACGAGGTCCACA
ACCN:	<u>BC049789, NM_001082476, NM_001082476.1, NM_001082476.2, BQ771207</u>
UniProt ID:	<u>A2AI05</u>
Synonyms:	4930447P04Rik; Ndor; NR1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)