

Product datasheet for **MP220268**

Nfyc Mouse qPCR Primer Pair (NM_001048168)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18046
Forward Sequence:	CCAGACACTTGCTACCAATGCC
Reverse Sequence:	AGGCATGGTGACTTGCTGGATC
ACCN:	<u>BC020117</u> , <u>BC053723</u> , <u>BC085261</u> , <u>NM_001048168</u> , <u>NM_001048168.1</u> , <u>NM_001048168.2</u> , <u>BE647422</u> , <u>BM461330</u> , <u>BU559204</u>
UniProt ID:	<u>P70353</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)