

Product datasheet for **MP220277**

Creb1 Mouse qPCR Primer Pair (NM_009952)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12912
Forward Sequence:	CACAGACCACTGATGGACAGCA
Reverse Sequence:	AGGACGCCATAACAACCTCCAGG
ACCN:	BC021649 , NM_009952 , NM_009952.1 , NM_009952.2 , BB831364 , BC048151 , BC058366 , BC071201 , BY315899 , BY782293
UniProt ID:	Q01147
Synonyms:	2310001E10Rik; 3526402H21Rik; AV083133; Creb; Creb-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)