

Product datasheet for **MP220295**

Slu7 Mouse qPCR Primer Pair (NM_198936)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	193116
Forward Sequence:	GAGTCTGTGAAGAAACCTCAGGC
Reverse Sequence:	TTCCGTTCTCATCGTCGCTGT
ACCN:	<u>NM_198936</u> , <u>NM_198936.1</u> , <u>BC013810</u> , <u>BC023057</u> , <u>BC025870</u> , <u>BC060954</u> , <u>BC082780</u> , <u>BC106099</u> , <u>BC148503</u> , <u>BY243619</u>
UniProt ID:	<u>Q8BHJ9</u>
Synonyms:	AU018913; D3Bwg0878e; D11Erttd730; D11Erttd730e
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)