

Product datasheet for **MP220380**

Sema3b Mouse qPCR Primer Pair (NM_009153)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20347
Forward Sequence:	CGAGCCTTCTTGGGACCTTTTG
Reverse Sequence:	GAAGGTGCCAAAGGTCTTGCTG
ACCN:	<u>BC090669</u> , <u>NM_009153</u> , <u>NM_009153.1</u> , <u>NM_009153.2</u> , <u>NM_009153.3</u> , <u>BB116052</u> , <u>BB787270</u> , <u>BC150758</u> , <u>BX515696</u> , <u>BY138852</u> , <u>BY341101</u>
Synonyms:	LUCA; Se; sem; SemA; sema5; Sema; semaV
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)