

Product datasheet for **MP220449**

Ang6 Mouse qPCR Primer Pair (NM_001011876)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	630952
Forward Sequence:	AAGCCAACAGGCTGGGATCACA
Reverse Sequence:	GCCTCCATTCTCTCCACAGATG
ACCN:	<u>NM_001011876</u> , <u>NM_001011876.1</u>
Synonyms:	EG630952
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)