

Product datasheet for **MP220496**

Mpl Mouse qPCR Primer Pair (NM_001122949)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	17480
Forward Sequence:	TCCTTGTAGAGGTGACCACAGC
Reverse Sequence:	ACTCCAGCCTTCCACTTGAGAC
ACCN:	NM_001122949 , NM_001122949.1 , NM_001122949.2 , BB850765 , BC101943 , BC103513 , BC103514 , BC103515 , BC146612 , BC146613
UniProt ID:	Q08351
Synonyms:	c-mpl; CD110; hlb219; TPO-R
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)