

Product datasheet for **MP220513**

Dab1 Mouse qPCR Primer Pair (NM_177259)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	13131
Forward Sequence:	GGCTCTATGTCTTTCGGCACTG
Reverse Sequence:	GGCGATCTGCTGTTGAACAAGG
ACCN:	NM_177259 , NM_177259.3 , NM_177259.4 , BC138538 , BC051993 , BC138539 , BE945099 , BE947086 , BI989333 , BY264959 , BY731247 , BY777583
UniProt ID:	P97318
Synonyms:	AI956902; C630028C02Rik; scm; scr; scrambler; yot
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)