

## Product datasheet for **MP220554**

### Ankle1 Mouse qPCR Primer Pair (NM\_172756)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	234396
Forward Sequence:	ACCAGAGGACACTATCTACGGC
Reverse Sequence:	TTCCTCATCCACCTGGCAGTCA
ACCN:	<u>NM_172756, NM_172756.1, NM_172756.2, NM_172756.3, BC030436, BC094316, BC145659</u>
UniProt ID:	<u>A8VU90</u>
Synonyms:	8430438L13Rik; Ankrd41; C430003I23; Lem3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)