

## Product datasheet for **MP220603**

### Ttc8 Mouse qPCR Primer Pair (NM\_198311)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	76260
Forward Sequence:	CAGAACTGCCTACACAGCACGA
Reverse Sequence:	CCTTTGCCAACTTAGGCTTCTGG
ACCN:	<u><a href="#">BC017523</a></u> , <u><a href="#">NM_198311</a></u> , <u><a href="#">NM_198311.1</a></u> , <u><a href="#">BY144447</a></u> , <u><a href="#">NM_198311.2</a></u>
UniProt ID:	<u><a href="#">Q8VD72</a></u>
Synonyms:	0610012F22Rik; AV001447; BBS8
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)