

Product datasheet for **MP220771**

Chd2 Mouse qPCR Primer Pair (NM_001081345)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	244059
Forward Sequence:	TCCGAAAGCAGGCACTGGATCA
Reverse Sequence:	GATGACGACTGTGCCGCTGAA
ACCN:	<u>NM_001081345, NM_001081345.1, NM_001081345.2, BC029703, BC062275</u>
UniProt ID:	<u>E9PZM4</u>
Synonyms:	2810013C04Rik; 2810040A01Rik; 5630401D06Rik; AI851092; BC029703; CHD-2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)