

Product datasheet for **MP220776**

Jak2 Mouse qPCR Primer Pair (NM_001048177)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	16452
Forward Sequence:	GCTACCAGATGGAAACTGTGCG
Reverse Sequence:	GCCTCTGTAATGTTGGTGAGATC
ACCN:	BC054807 , BC059834 , NM_001048177 , NM_001048177.1 , NM_001048177.2 , BE981299 , BF019856 , BF730796 , BG068603
UniProt ID:	Q62120
Synonyms:	Fd17
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)