

Product datasheet for **MP220951**

Gcnt2 Mouse qPCR Primer Pair (NM_133219)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	14538
Forward Sequence:	CATTGCTCTCTTAGAGTGGTCTAA
Reverse Sequence:	CCATCCACTTCACAGCTCTGAG
ACCN:	BC054845 , BC094572 , NM_133219 , NM_133219.1 , BE989476 , BM933481 , BQ892107 , BY068011 , BY208177
UniProt ID:	P97402
Synonyms:	5330430K10Rik; IGnT; IGnTA; IGnTB; IGnTC
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)