

## Product datasheet for **MP220961**

### **Txn14a Mouse qPCR Primer Pair (NM\_025299)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	27366
<b>Forward Sequence:</b>	GTGGTCGTCATTCGTTTCGGAC
<b>Reverse Sequence:</b>	CGTACATCTTGTTGAAGTCAGGC
<b>ACCN:</b>	<a href="#">BC014770</a> , <a href="#">BC031634</a> , <a href="#">NM_025299</a> , <a href="#">NM_025299.1</a> , <a href="#">NM_025299.2</a> , <a href="#">NM_025299.3</a> , <a href="#">BC003876</a> , <a href="#">BC022582</a> , <a href="#">BC025053</a> , <a href="#">BC029327</a> , <a href="#">BC057058</a> , <a href="#">BG074208</a> , <a href="#">BG086914</a> , <a href="#">BU514237</a> , <a href="#">BY065764</a> , <a href="#">BY098388</a> , <a href="#">BY654643</a> , <a href="#">NM_025299.4</a>
<b>UniProt ID:</b>	<a href="#">P83877</a>
<b>Synonyms:</b>	D18Wsu98e; Dim1; ENSMUSG00000057130; Txn14; U5-15kD; U5-15kDa
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)