

Product datasheet for **MP220990**

Adora1 Mouse qPCR Primer Pair (NM_001039510)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11539
Forward Sequence:	GATCGGTACCTCCGAGTCAAGA
Reverse Sequence:	CACTCAGGTTGTTCCAGCCAAAC
ACCN:	<u>BC079624</u> , <u>NM_001039510</u> , <u>NM_001039510.1</u> , <u>NM_001039510.2</u> , <u>BC026602</u> , <u>BC058421</u> , <u>BY257408</u> , <u>BY737264</u>
UniProt ID:	<u>Q60612</u>
Synonyms:	A1-AR; A1AR; A1R; AA1R; AI848715; BB176431; Ri
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)