

## Product datasheet for **MP220998**

### Fiz1 Mouse qPCR Primer Pair (NM\_001110328)

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | qPCR Primer Pairs  |
| Gene ID:          | 23877  |
| Forward Sequence: | GCTTCAACTTGGCTAACCACCTG  |
| Reverse Sequence: | GGCTTTTCACCAGTGTGGACCA   |
| ACCN:             | <a href="#">BC006633</a> , <a href="#">BC106113</a> , <a href="#">NM_001110328</a> , <a href="#">NM_001110328.1</a> , <a href="#">BB855369</a> , <a href="#">BC138231</a> , <a href="#">BC138232</a> , <a href="#">BC145207</a> , <a href="#">BU053482</a> , <a href="#">BY344494</a> , <a href="#">NM_001110328.2</a>                       |
| UniProt ID:       | <a href="#">Q9WTJ4</a>   |
| Synonyms:         | AI790204; mFLJ00416  |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.  |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.  |
| Stability:        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)