

Product datasheet for **MP221035**

Rps6kb1 Mouse qPCR Primer Pair (NM_001114334)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	72508
Forward Sequence:	AGGTGGAACCTCCCTTAAGCC
Reverse Sequence:	CCAGAAAGACCTGGTTGGCACT
ACCN:	<u>BC038491</u> , <u>BC071244</u> , <u>NM_001114334</u> , <u>NM_001114334.1</u>
UniProt ID:	<u>Q8BSK8</u>
Synonyms:	70kDa; 2610318I15Rik; 4732464A07Rik; AA959758; AI256796; AI314060; p70/85s6k; p70s6k; S6K1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)