

Product datasheet for **MP221092**

Nol4l Mouse qPCR Primer Pair (NM_001134300)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	329540
Forward Sequence:	CATCCGCACATACCTCAAGTCC
Reverse Sequence:	GCATCCTCTTGGCTGCCTTTCT
ACCN:	<u>BC075648</u> , <u>NM_001134300</u> , <u>NM_001134300.1</u> , <u>NM_001134300.2</u> , <u>BB850016</u> , <u>BC029060</u> , <u>BC079552</u> , <u>BC094234</u>
Synonyms:	4931410H15; 8430427H17Rik; AI429157; Gm1636
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)